



ENDURE MEDICAL, INC.


Suite 1601, 16th Floor Tekliffe East Tower, Exchange Road, Ortigas Center
Brgy. San Antonio, Pasig City, Philippines, 1605

ENDURE Surgical Facemask, Earloop



INTENDED USE

Intended to be used to protect the patient from infective agents by covering the mouth and nose of medical staff. It may also be intended to be worn by patients and other persons to reduce the risk of spread of infections, particularly in epidemic or pandemic situations.



SPECIFICATION

- **Model No. :** N77A02
- **Material:** PP + Filter Paper + PP
 - Made of non-woven and melt-blown fabric
- **Weight:** 25gsm + 25 filter paper + 25gsm
- 3-ply
- With flexible nose wire
- Latex free and fiberglass free.
- **Size:** 17.5 cm x 9.5 cm
- Non-sterile
- Disposable
- Bacterial Filtration Efficiency $\geq 98\%$, Type IIR
- **Color:** Blue and White
- **Packaging:** 50pcs/bag/box, 2000pcs/ctn



FDA CMDN No.: CMDN-2025-1126758



SHELF LIFE: 3 Years



STORAGE CONDITION:

- (1) Store in dry, clean conditions in original packaging.
- (2) Store away from direct sunlight, source of high temperature and solvent vapors.
- (3) Store with the temperature range -5°C to $+45^{\circ}\text{C}$ and with relative humidity below 65%.